

HEDY Swan Digital Mammography

Items	Configuration 1	Configuration 2 - <i>Tomosynthesis</i>
Tube	IAE tube 300KHU	IAE tube 300KHU
Detector	ANRAD amorphous selenium flat panel 24*30 cm	HEDY 2430FD amorphous selenium flat panel 24*30 cm
Generator	Spellman	Spellman
Rack	HEDY C-arm rack	HEDY C-arm rack
Multifunctional workstation	HEDY multifunctional breast DR system workstation	HEDY multifunctional breast DR system workstation